



Antigenic integral membrane glycoprotein Antibody

Product Code	CSB-PA328381ZA01SXQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23126
Immunogen	Recombinant Schistosoma mansoni Antigenic integral membrane glycoprotein protein
Raised In	Rabbit
Species Reactivity	Schistosoma mansoni
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Schistosoma mansoni (Blood fluke)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Antigenic integral membrane glycoprotein,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.