



pht3 Antibody

Product Code	CSB-PA313137ZA01FFZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05183
Immunogen	Recombinant Pseudomonas putida pht3 protein
Raised In	Rabbit
Species Reactivity	Pseudomonas putida
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pseudomonas putida (Arthrobacter siderocapsulatus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pht3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.