



yafQ Antibody

Product Code	CSB-PA332059ZA01HTA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P44041
Immunogen	Recombinant Haemophilus influenzae yafQ protein
Raised In	Rabbit
Species Reactivity	Haemophilus influenzae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
Lead Time	Made-to-order (12-14 weeks)
Target Names	yafQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.