



E & M Antibody

Product Code	CSB-PA34028ZA01GND
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DTC4&P0DTC5
Immunogen	Recombinant Severe acute respiratory syndrome coronavirus 2 E & M protein
Raised In	Rabbit
Species Reactivity	Severe acute respiratory syndrome coronavirus 2
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	E & M