



H3L Antibody

Product Code	CSB-PA324943ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20497
Immunogen	Recombinant Vaccinia virus H3L protein
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	H3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.