



HA Antibody

Product Code	CSB-PA321868ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20978
Immunogen	Recombinant Vaccinia virus HA protein
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	HA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.