



ibpA Antibody

Product Code	CSB-PA619558ZA01YAG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1CD74
Immunogen	Recombinant Yersinia pestis bv. Antiqua ibpA protein
Raised In	Rabbit
Species Reactivity	Yersinia pestis bv. Antiqua
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Yersinia pestis bv. Antiqua (strain Nepal516)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ibpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.