



cspA Antibody

Product Code	CSB-PA359368ZA01SBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A9Y5
Immunogen	Recombinant Salmonella enteritidis cspA protein
Raised In	Rabbit
Species Reactivity	Salmonella enteritidis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Salmonella enteritidis
Lead Time	Made-to-order (12-14 weeks)
Target Names	cspA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.