



P2 Antibody

Product Code	CSB-PA319967ZA01PUV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11124
Immunogen	Recombinant Pseudomonas phage phi6 P2 protein
Raised In	Rabbit
Species Reactivity	Pseudomonas phage phi6
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pseudomonas phage phi6 (Bacteriophage phi-6)
Lead Time	Made-to-order (12-14 weeks)
Target Names	P2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.