



cla Antibody

Product Code	CSB-PA316109ZA01AVC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CL71
Immunogen	Recombinant Aspergillus clavatus cla protein
Raised In	Rabbit
Species Reactivity	Aspergillus clavatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Aspergillus clavatus (strain ATCC 1007 / CBS 513.65 / DSM 816 / NCTC 3887 / NRRL 1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	cla
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.