



SAP2 Antibody

Product Code	CSB-PA505575ZA01CZE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C4YMJ3
Immunogen	Recombinant Candida albicans SAP2 protein
Raised In	Rabbit
Species Reactivity	Candida albicans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Candida albicans (strain WO-1) (Yeast)
Lead Time	Made-to-order (12-14 weeks)
Target Names	SAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.