



Major capsid protein VP1 Antibody

Product Code	CSB-PA355947ZA01HBZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03092
Immunogen	Recombinant Hamster polyomavirus Major capsid protein VP1 protein
Raised In	Rabbit
Species Reactivity	Hamster polyomavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Hamster polyomavirus (HaPyV) (Mesocricetus auratus polyomavirus 1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Major capsid protein VP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.