



PRKACA Antibody

Product Code	CSB-PA822648ZA01FQF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8MJ44
Immunogen	Recombinant Canis lupus familiaris PRKACA protein
Raised In	Rabbit
Species Reactivity	Canis lupus familiaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PRKACA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.