



Early E3 18.5 kDa glycoprotein Antibody

Product Code	CSB-PA301615ZA01HIK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68978
Immunogen	Recombinant Human adenovirus C serotype 2 Early E3 18.5 kDa glycoprotein protein
Raised In	Rabbit
Species Reactivity	Human adenovirus C serotype 2
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human adenovirus C serotype 2 (HAdV-2) (Human adenovirus 2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Early E3 18.5 kDa glycoprotein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.