



Dihydrolipoyl dehydrogenase Antibody

Product Code	CSB-PA006928ZA01HRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80647
Immunogen	Recombinant Hymenolepis diminuta Dihydrolipoyl dehydrogenase protein
Raised In	Rabbit
Species Reactivity	Hymenolepis diminuta
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Hymenolepis diminuta (Rat tapeworm)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Dihydrolipoyl dehydrogenase
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.