



DODA Antibody

Product Code	CSB-PA310823ZA01AZZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87064
Immunogen	Recombinant Amanita muscaria DODA protein
Raised In	Rabbit
Species Reactivity	Amanita muscaria
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Amanita muscaria (Fly agaric)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	DODA