



7b Antibody

Product Code	CSB-PA756006ZA01HQE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7TFA1
Immunogen	Recombinant Human SARS coronavirus 7b protein
Raised In	Rabbit
Species Reactivity	Human SARS coronavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human SARS coronavirus (SARS-CoV) (Severe acute respiratory syndrome coronavirus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	7b
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.