




cry1Ia Antibody

Product Code	CSB-PA673817ZA01BDC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q45752
Immunogen	Recombinant Bacillus thuringiensis subsp. Kurstaki Pesticidal crystal protein cry1Ia(cry1Ia),partial (97-180aa)
Species Reactivity	Bacillus thuringiensis subsp. Kurstaki
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus thuringiensis subsp. Kurstaki
Target Names	cry1Ia

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: cry1Ia antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36.2 kDa

Observed band size: 37 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.