



CAT Antibody

Product Code	CSB-PA881920ZA01DQD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9M5L6
Immunogen	Recombinant Capsicum annuum CAT protein
Raised In	Rabbit
Species Reactivity	Capsicum annuum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Capsicum annuum (Bell pepper)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	CAT