



MTR_4g091010 Antibody

Product Code	CSB-PA34528ZA01MQP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G7JRT5
Immunogen	Recombinant Medicago truncatula MTR_4g091010 protein
Raised In	Rabbit
Species Reactivity	Medicago truncatula
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Medicago truncatula (Barrel medic) (Medicago tribuloides)
Lead Time	Made-to-order (12-14 weeks)
Target Names	MTR_4g091010
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.