



Genome polyprotein Antibody

Product Code	CSB-PA22398ZA01GLD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A0N9XAY5
Immunogen	Recombinant Genome polyprotein protein
Raised In	Rabbit
Species Reactivity	dengue virus type 2
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	dengue virus type 2
Lead Time	Made-to-order (12-14 weeks)
Target Names	Genome polyprotein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.