



clp1 Antibody

Product Code	CSB-PA740650ZA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NS21
Immunogen	Recombinant Xenopus laevis clp1 protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Xenopus laevis (African clawed frog)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	clp1