



nbIA Antibody

Product Code	CSB-PA333966ZA01FPY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35087
Immunogen	Recombinant Synechococcus elongatus nbIA protein
Raised In	Rabbit
Species Reactivity	Synechococcus elongatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Synechococcus elongatus (strain PCC 7942 / FACHB-805) (Anacystis nidulans R2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	nbIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.