



pucL Antibody

Product Code	CSB-PA523373ZA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O32141
Immunogen	Recombinant Bacillus subtilis pucL protein
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pucL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.