



psmA1 Antibody

Product Code	CSB-PA314123ZA01FLF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C7Y1
Immunogen	Recombinant Staphylococcus aureus psmA1 protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Staphylococcus aureus (strain NCTC 8325 / PS 47)
Lead Time	Made-to-order (12-14 weeks)
Target Names	psmA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.