



JBP1 Antibody

Product Code	CSB-PA390102ZA01LOF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4HU70
Immunogen	Recombinant Leishmania infantum JBP1 protein
Raised In	Rabbit
Species Reactivity	Leishmania infantum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Leishmania infantum
Lead Time	Made-to-order (12-14 weeks)
Gene Names	JBP1