



Capsid protein VP2 Antibody

Product Code	CSB-PA351680ZA01CCR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61826
Immunogen	Recombinant Canine parvovirus type 2 Capsid protein VP2 protein
Raised In	Rabbit
Species Reactivity	Canine parvovirus type 2
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Canine parvovirus type 2 (CPV-2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Capsid protein VP2,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.