



FKTN Antibody

Product Code	CSB-PA008709ZA01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q60HG0
Immunogen	Recombinant Macaca fascicularis FKTN protein
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	FKTN