



SIGLEC1 Antibody

Product Code	CSB-PA021293ZA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7LCJ3
Immunogen	Recombinant Sus scrofa SIGLEC1 protein
Raised In	Rabbit
Species Reactivity	Sus scrofa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Sus scrofa (Pig)
Lead Time	Made-to-order (12-14 weeks)
Target Names	SIGLEC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.