



cyp Antibody

Product Code	CSB-PA35408ZA01GME
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q50927
Immunogen	Recombinant Nitrosomonas europaea cyp protein
Raised In	Rabbit
Species Reactivity	Nitrosomonas europaea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Nitrosomonas europaea (bioreactor metagenome)
Lead Time	Made-to-order (12-14 weeks)
Target Names	cyp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.