



IL7 Antibody

Product Code	CSB-PA011669ZA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26895
Immunogen	Recombinant Bos taurus IL7 protein
Raised In	Rabbit
Species Reactivity	Bos taurus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bos taurus (Bovine)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IL7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.