



5b Antibody

Product Code	CSB-PA323968ZA01ARV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19745
Immunogen	Recombinant Avian infectious bronchitis virus 5b protein
Raised In	Rabbit
Species Reactivity	Avian infectious bronchitis virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Avian infectious bronchitis virus (strain KB8523) (IBV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	5b