



inIA Antibody

Product Code	CSB-PA338647ZA01LPY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DJM0
Immunogen	Recombinant Listeria monocytogenes serovar 1/2a inIA protein
Raised In	Rabbit
Species Reactivity	Listeria monocytogenes serovar 1/2a
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Listeria monocytogenes serovar 1/2a (strain ATCC BAA-679 / EGD-e)
Lead Time	Made-to-order (12-14 weeks)
Target Names	inIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.