



ftsQ Antibody

Product Code	CSB-PA352568ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WNA0
Immunogen	Recombinant Mycobacterium tuberculosis ftsQ protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium tuberculosis
Lead Time	Made-to-order (12-14 weeks)
Gene Names	ftsQ