



proC Antibody

Product Code	CSB-PA607537ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WHU7
Immunogen	Recombinant Mycobacterium tuberculosis proC protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium tuberculosis
Lead Time	Made-to-order (12-14 weeks)
Target Names	proC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.