



sodC Antibody

Product Code	CSB-PA022397ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WGE8
Immunogen	Recombinant Mycobacterium tuberculosis sodC protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium tuberculosis
Lead Time	Made-to-order (12-14 weeks)
Target Names	sodC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.