



CCL17 Antibody

Product Code	CSB-PA856825ZA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q95N01
Immunogen	Recombinant Canis familiaris CCL17 protein
Raised In	Rabbit
Species Reactivity	Canis familiaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Canis familiaris (Dog) (Canis lupus familiaris)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CCL17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.