



Genome polyprotein Antibody

Product Code	CSB-PA36438ZA01GOZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q32ZE1
Immunogen	Recombinant Zika virus Genome polyprotein protein
Raised In	Rabbit
Species Reactivity	Zika virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Zika virus (ZIKV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Genome polyprotein,partial