



Lectin Antibody

Product Code	CSB-PA355890ZA01LLM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02870
Immunogen	Recombinant Lens culinaris Lectin protein
Raised In	Rabbit
Species Reactivity	Lens culinaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Lens culinaris (Lentil) (Cicer lens)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Lectin
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.