



bipD Antibody

Product Code	CSB-PA663187ZA01BAAL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3JL26
Immunogen	Recombinant Burkholderia pseudomallei bipD protein
Raised In	Rabbit
Species Reactivity	Burkholderia pseudomallei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Burkholderia pseudomallei (strain 1710b)
Lead Time	Made-to-order (12-14 weeks)
Target Names	bipD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.