



aiiA Antibody

Product Code	CSB-PA315623ZA01BDC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CJ63
Immunogen	Recombinant Bacillus thuringiensis subsp. Kurstaki aiiA protein
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. Kurstaki
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bacillus thuringiensis subsp. Kurstaki
Lead Time	Made-to-order (12-14 weeks)
Gene Names	aiiA