



kgp Antibody

Product Code	CSB-PA690409ZA01PQP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q51817
Immunogen	Recombinant Porphyromonas gingivalis kgp protein
Raised In	Rabbit
Species Reactivity	Porphyromonas gingivalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Porphyromonas gingivalis
Lead Time	Made-to-order (12-14 weeks)
Target Names	kgp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.