



PHLPVI Antibody

Product Code	CSB-PA337347ZA01EUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43215
Immunogen	Recombinant Phleum pratense PHLPVI protein
Raised In	Rabbit
Species Reactivity	Phleum pratense
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Phleum pratense (Common timothy)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	PHLPVI