



ideZ Antibody

Product Code	CSB-PA35978ZA01GOM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0PIW1
Immunogen	Recombinant Streptococcus equi subsp. zooepidemicus ideZ protein
Raised In	Rabbit
Species Reactivity	Streptococcus equi subsp. zooepidemicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Streptococcus equi subsp. zooepidemicus
Lead Time	Made-to-order (12-14 weeks)
Target Names	ideZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.