



5a Antibody

Product Code	CSB-PA848658ZA01BJF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8V434
Immunogen	Recombinant Bovine coronavirus 5a protein
Raised In	Rabbit
Species Reactivity	Bovine coronavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bovine coronavirus (strain 98TXSF-110-LUN) (BCoV-LUN) (BCV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	5a