



Goat Anti-Human IgG, Fc γ fragment specific; FITC conjugated

Product Code	CSB-PA271513
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Human
Tested Applications	FC; Recommended dilution: FC:1:50-1:200
Form	Liquid
Conjugate	FITC
Storage Buffer	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Purification Method	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Isotype	IgG
Clonality	Polyclonal
Product Type	Secondary Antibody
Immunogen Species	Homo sapiens (Human)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.