



Goat Anti-Mouse IgG, Fc γ fragment specific; HRP conjugated

Product Code	CSB-PA209995
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Mouse
Tested Applications	ELISA, WB, IHC, ICC; Recommended dilution: ELISA:1:500-1:10000, WB:1:3000-1:5000, IHC:1:50-1:500, ICC:1:50-1:500
Form	Liquid
Conjugate	HRP
Storage Buffer	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Purification Method	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Isotype	Secondary Antibody
Product Type	Secondary Antibody
Immunogen Species	Mouse