




HEK293F/Human DLL3 Stable Cell Line

Product Code	CSB-SC882142HU
Target Name	suspension
Storage	Frozen in liquid nitrogen or stored at -80°C
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	90% HEK293-F Expression Medium+10% DMSO
Resistance	Puromycin(2 µg/mL)
Intended Use	HEK293-F Expression Medium

Image



Untransfected HEK293 cells (green line) and transfected Human DLL3 HEK293 stable cells (red line) were stained with anti-DLL3 antibody (rovalpituzumab-like) (2µg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Human IgG1 Fc antibody and analyzed with flow cytometry.