




HEK293T/Human SEMA4D Stable Cell Line

Product Code	CSB-SC835707HU
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	90% FBS+10% DMSO
Resistance	Puromycin(2 µg/mL)
Intended Use	DMEM medium +10% FBS

Image



Untransfected HEK293T cells (green line) and transfected Human SEMA4D HEK293T stable cells (red line) were stained with anti-SEMA4D antibody (CSB-RA835707A2HU) (2µg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.