




African swine fever virus Phosphoprotein p30(Ba71V-93) Monoclonal Antibody

Product Code	CSB-MA335904A1m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34204
Immunogen	Recombinant Phosphoprotein p30 protein (1-204AA)
Raised In	Mouse
Species Reactivity	African swine fever virus
Tested Applications	ELISA
Relevance	Modifies the subcellular distribution of heterogeneous nuclear ribonucleoprotein K (HNRNPK) and may contribute to modulate HNRNPK functions related to processing and export of mRNAs during ASFV infection.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG1
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Immunogen Species	African swine fever virus
Research Area	Others
Target Names	Ba71V-93
Clone No.	5F2F3

Image



The Binding Activity of African swine fever virus Phosphoprotein p30(Ba71V-93). Activity: Measured by its binding ability in a functional ELISA. Immobilized African swine fever virus Phosphoprotein p30(Ba71V-93) protein at 2 µg/ml can bind African swine fever virus Phosphoprotein p30(Ba71V-93) Antibody, the EC₅₀ is 0.01949 to 0.02692 ug/ml.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.